



# Jazz Sensor Software Compatibility Chart

Dental Imaging Software	Direct Integration	TWAIN
ACTEON	X	✓
AEKA	✓	✓
ARCHY	X	✓
CARESTREAM	✓	✓
CDR DICOM	✓	✓
CLIOSOFT	X	✓
CURVE	X	✓
DENTTIO	X	X
DENTRIX ASCEND	✓	✓
DENTRIX IMAGING	✓	✓
DEXIS 9	✓	✓
DEXIS 10	✓	✓

DEXIS 11 - V11.2.2.18 and earlier	✓	✓
DEXIS 11 - after V11.2.2.18	✗	✓
DTX STUDIO	✗	✓
EAGLESOFT	✓	✓
EZDENT-i	✗	✓
GX PICTURE V3.0.0 to V3.8.0	✓	✓
GX PICTURE V4.1, V4.2	✗	✗
iDENTALSOF	✗	✓
MIPACS	✗	✓
MOGO	✗	✗
ORYX	✓	✓
OPEN DENTAL	✓	✓
OVERJET IRIS	✓	✓
PATIENT GALLERY	✗	✓
ROMEXIS	✗	✓

SIDEXIS V4.1.x - 4.2.x	✓	✓
SIDEXIS 4.3.x and above	X	✓
SOTA	X	✓
TAB32	X	✓
TDO	X	✓
TIGERVIEW	X	✓
TRACKER	X	✓
VISIX	X	✓
VIXWIN	✓	✓
XDR	X	✓
XV WEB	✓	✓