



Jazz Sensor Software Compatibility Chart

Dental Imaging Software	Direct Integration	TWAIN
ACTEON	X	✓
CARESTREAM	✓	✓
CDR DICOM	✓	✓
CLIOSOFT	X	✓
CURVE	X	✓
DENTTIO	X	X
DENTRIX ASCEND	✓	✓
DENTRIX IMAGING	✓	✓
DEXIS 9	✓	✓
DEXIS 10	✓	✓
DEXIS 11 - V11.2.2.18 and earlier	✓	✓
DEXIS 11 - after V11.2.2.18	X	✓

DTX STUDIO	X	✓
EAGLESOFT	✓	✓
GX PICTURE V3.0.0 to V3.8.0	✓	✓
GX PICTURE V4.1, V4.2	X	X
MIPACS	X	✓
MOGO	X	X
ORYX	✓	✓
OPEN DENTAL	✓	✓
PATIENT GALLERY	X	✓
ROMEXIS	X	✓
SIDEXIS	✓	✓
SOTA	X	✓
TAB32	X	✓
TIGERVIEW	X	✓
TRACKER	X	✓

VISIX	X	✓
VIXWIN	✓	✓
XDR	X	✓
XV WEB	✓	✓