

Jazz Sensor Software Compatibility Chart

Dental Imaging Software	Direct Integration	TWAIN
ACTEON	X	✓
AEKA	✓	✓
APTERYX	✓	✓
ARCHY	X	✓
CARESTREAM	✓	✓
CDR DICOM	✓	✓
CLIOSOFT	X	✓
CURVE	X	✓
DENTTIO	X	X
DENTRIX ASCEND	✓	✓
DENTRIX IMAGING	✓	✓
DEXIS 9	✓	✓

DEXIS 10	✓	✓
DEXIS 11 - V11.2.2.18 and earlier	✓	✓
DEXIS 11 - after V11.2.2.18	X	✓
DTX STUDIO	X	✓
EAGLESOFT	✓	✓
EZDENT-i	X	✓
GX PICTURE	✓	✓
iDENTALSOF	X	✓
MIPACS	✓	✓
MOGO	X	X
ORYX	✓	✓
OPEN DENTAL (DESKTOP)	X	✓
OPEN DENTAL (CLOUD)	X	✓
OVERJET IRIS	✓	✓

PATIENT GALLERY	X	✓
ROMEXIS	X	✓
SIDEXIS	✓	✓
SOTA CLOUD	✓	✓
TAB32	X	✓
TDO	X	✓
TIGERVIEW	X	✓
TRACKER	X	✓
VISIX	X	✓
VIXWIN	✓	✓
XDR	X	✓